

Batteries: What You Need to Know

NOLA HEARING



Battery sizes and colors are universal. This means that no matter where you are when you run out of batteries, you can purchase new ones as long as you remember either the color or the number.

The corresponding numbers and colors are as follows:



10, yellow



13, orange



312, brown



675, blue

There are a few exceptions to this rule. Occasionally, there will be a smaller or larger lip on the batteries that will prevent the contacts inside your hearing aid from meeting, thus keeping the instrument from working properly or at all.

If you find that you have purchased the incorrect size batteries or you are having another issue, please bring them to us and we will exchange them for you at no additional charge.